

RISVEGLIO CUTANEO: UN BOOSTER DI NUOVA GENERAZIONE IN PROTOCOLLO COMBINATO

Elena Torre

OBIETTIVI

Valutazione dell'effetto Anti-aging con un nuovo protocollo bio-ristrutturante e bio-rivitalizzante, nato dalla combinazione di un trattamento iniettivo con un nuovo Booster a base di HA debolmente cross-linkato e un soft-peeling con un mix di 12 attivi sinergici per il miglioramento della skin quality.

MATERIALI E METODI

IL NOSTRO SKIN BOOSTER

Il Booster utilizzato, è un trattamento iniettivo a base di acido ialuronico monofasico a bassa reticolazione, indicato per cute asfittica/disidratata, correzione delle micro-rugosità superficiali, miglioramento dell'idratazione e turgore cutaneo.

Utilizzato in protocollo integrato con un soft-peeling di nuova generazione contenente un mix di 12 attivi sinergici, a integrare l'azione profonda dello skinbooster.

SKIN BOOSTER RIASSORBIBILE che si diffonde con naturalezza. Viscosità ed elasticità dell'Acido Ialuronico sono tali da conservare la coesione del volume senza però ostacolarne la capacità di adattarsi alla morfologia individuale.

- Tecnica a microboli (5 punti bioestetici di trazione)
- Tecnica a ventaglio (ago o microcannula) per gli inestetismi più evidenti

SOFT-PEELING IN SISTEMA INTEGRATO IN DUE FASI PEELING indicato per riattivare il turn over epidermico, rendere più uniforme la pigmentazione e aumentare la luminosità. Si compone di un peeling in quattro step: due fasi peeling con differenti acidi a diverse concentrazioni, due fasi di recovery

PROCEDURA

PROTOCOLLO COMBINATO NELLA STESSA SEDUTA

- SOFT-PEELING + SKIN BOOSTER: 1 ogni 3 settimane per 2-3 sedute

PROTOCOLLO ALTERNATO

Alternare i due trattamenti singoli ogni 15gg, per 2-3 sedute:

- SOFT-PEELING
- SKIN BOOSTER

RISULTATI

Il protocollo è risultato essere molto efficace nella riduzione dei danni da aging, come bio-ristrutturante, bio-rivitalizzante, con leggero effetto liftante.

La combinazione ha contribuito ad attenuare le discromie, aumentare la luminosità e nelle pelli seborroiche, svolgere un'azione astringente e purificante.

CONCLUSIONI

Il protocollo è risultato efficace nella riduzione dei danni da foto e crono aging, con leggero effetto liftante. Versatile, con un'alta compliance e soddisfazione da parte dei pazienti

SKIN AWAKENING: A NEW GENERATION BOOSTER IN A COMBINED PROTOCOL

OBJECTIVES

Evaluation of the Anti-aging effect with a new bio-restructuring and bio-revitalizing protocol, born from the combination of an injection treatment with a new Booster based on HA weakly cross-linked and a soft-peeling with a mix of 12 synergic active ingredients for the improvement of skin quality.

MATERIALS AND METHODS

OUR SKIN BOOSTER

The Booster used is an injection treatment based on monophasic hyaluronic acid with low cross-linking, indicated for dehydrated skin, correction of superficial micro-rugosity, improvement of hydration and skin turgidity.

Used in integrated protocol with a new generation soft-peeling containing a mix of 12 synergic active ingredients, to integrate the deep action of the skinbooster.

Re-absorbable SKIN BOOSTER that spreads naturally. Viscosity and elasticity of the Hyaluronic Acid are such as to preserve the cohesion of the volume without hindering its ability to adapt to individual morphology.

- Microbolus technique (5 bio-aesthetic traction points)
- Fan technique (needle or microcannula) for the most evident imperfections

SOFT-PEELING IN TWO-PHASE INTEGRATED SYSTEM PEELING indicated to reactivate the epidermal turn over, make pigmentation more uniform and increase the brightness. It consists of a four-step peeling: two phases of peeling with different acids at different concentrations, two phases of recovery, and one phase of recovery.

PROCEDURE

COMBINED PROTOCOL IN THE SAME SESSION

- SOFT-PEELING + SKIN BOOSTER: 1 every 3 weeks for 2-3 sessions

ALTERNATE PROTOCOL

Alternate the two single treatments every 15 days for 2-3 sessions:

- SOFT-PEELING
- SKIN BOOSTER

RESULTS

The protocol was found to be very effective in reducing aging damage, as a bio-restructuring, bio-revitalizing, with a slight lifting effect.

The combination has contributed to attenuate dyschromia, increase the brightness and in seborrhoeic skins, carry out an astringent and purifying action.

CONCLUSIONS

The protocol was effective in the reduction of photo and chrono aging damages, with a slight lifting effect. Versatile, with high compliance and satisfaction from patients.

1. Syed Nasir Abbas Bukhari et al. Int J Biol Macromol. 2018 Dec.

Hyaluronic acid, a promising skin rejuvenating biomedicine: A review of recent updates and pre-clinical and clinical investigations on cosmetic and nutricosmetic effects

2. Conciliun Feel IFU – Conciliun Feel Booster

3. Data on file - Valutazione clinica e strumentale dell'efficacia del trattamento con B-TWIN SYSTEM per il miglioramento della skin-quality